promed EMT-4

Promed GmbH

Lindenweg 11

D-82490 Farchant

Tel: +49 (0)8821/9621-0

Fax: +49 (0)8821/9621-21

info@promed.de

www.promed.de

Gebrauchsanleitung | Instruction leaflet | Instruction d'utilisation | Istruzioni per l'uso | Instrucciones de empleo

promed EMT-4

Digitales, elektrisches Schmerztherapie- und Muskelstimulations-Gerät (TENS & EMS)
(TENS: Transkutane Elektrische Nerven-Stimulation) (EMS: Elektrische Muskel-Stimulation)

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








PRODUCT INFORMATION / ACCESSORIES


About this product:

Your **Promed EMT-4** is a battery-operated unit with two adjustable output channels. This combined TENS-EMS unit generates electrical pulses whose amplitude, duration and modulation can be varied via the buttons. The **Promed EMT-4** can be prescribed by your doctor (classification number 09.31.02.1-999).

- 1** Output sockets for electrode cables
- 2** On/Off button
- 3** Program select button
- 4** Intensity regulation buttons
- 5** Confirm button
- 6** Channel select button
- 7** Battery compartment with 4 x 1,5 V batteries, AAA type
- 8** Electrode plugs (cable – electrodes)
- 9** Connector plugs (cable – main unit)

Product features:


- 1** Eleven preset programs (see program table) and two programmable programs
- 2** LCD (Liquid Crystal Display) for: **a)** Treatment time; **b)** Program; **c)** Frequency (pulses per second); **d)** Pulse width (duration of individual pulses); **e)** Intensity
- 3** Two channels, independently adjustable in 50 stages
- 4** Connection check: when the electrodes are not connected or attached, the intensity cannot be set higher than stage 1, or reverts to "0". (to avoid draining the batteries)
- 5** Application protection: when the unit is switched on, the intensity is set to stage 15 (of 50).
- 6** Memory function: the unit stores the treatment duration.
- 7** Time display: displays the time in "sleep mode".



PACKAGE CONTENTS

System components

Your **Promed EMT-4**-unit includes the following accessories::

- ❶ TENS/EMS main unit: 1
- ❷ Electrode connector cables: 2
- ❸ Self-adhesive electrodes: 2 pairs (4 pads)
- ❹ Unit stand: 1
- ❺ Operating instructions: 1
- ❻ Batteries: 4
- ❼ Storage case: 1

INTRODUCTION

What is pain?

Pain is the body's own early warning system. Feeling pain is important, as it indicates an abnormal condition within the body and serves as a warning to us before any further damage or injury occurs. However, long-term, persistent pains – often known as chronic pain – serve no obvious useful purpose once a diagnosis has been made. TENS was developed to soothe or eliminate certain types of chronic and acute pain.

We differentiate between two fundamental types of pain:


- **Acute pain**

as a chief symptom can often help the doctor with diagnosis and the acute pain has a protective function for the patient.

- **Chronic pain**

can often become part of the illness in itself. A patient suffering from chronic pain will often suffer for years and experiences changes in his/her personality structure.

The term "stimulation current therapy" serves as an umbrella term for the electrical stimulation of nerves and muscles.

Basic terms in electrical pain therapy

Pulse frequency

Number of pulses per second,
measured in Hertz (Hz).
Tens works between 1 and 200 Hz.

Pulse width




The duration of an individual pulse,
measured in microseconds. 1µsec.
corresponds to 1/1.000.000 second.

Current intensity

Measured in milliamperes 1mA = 1/1000 A.
The mains outlet outputs 16 A.
Tens units work with a maximum of 50 mA.

Modulation

Pulse intensity, pulse frequency and pulse duration
are simultaneously varied.

INTRODUCTION

How does TENS work?

TENS is a non-invasive and drug-free method of treating pain. The TENS unit transmits soothing electrical pulses through the skin that stimulate the nerve (or nerves) in the area being treated. In most cases, this stimulation greatly reduces or even completely eliminates the sensation of pain, as it blocks the original pain signals to the brain.

It is also reported that TENS stimulation promotes the release of endorphins into the bloodstream, leading to a further reduction in pain.

The benefits of TENS devices in pain management have been clinically proven in the case of many patients. To get the most out of your TENS unit, read these instructions carefully and follow your doctor's recommendations for treatment.

All information, general sensations, skin sensations, deep sensibility, is relayed to the CNS (central nervous system) via the nerve pathways.

The classic TENS method **interrupts pain transmission** for an indefinite period of time.

It does not affect the source of pain locally, but acts **exclusively on the nerve tissue**.

Nerve tissue can be stimulated by electrical current **over a wide spectrum of frequencies**.

Because of the interaction between the stimulation frequency and the transmission speed of the nerve fibers, it can be seen **that different frequencies selectively activate certain types of nerve**.

The body has two different types of nerve fiber:

- a) Unmyelinated nerve fibers with continuous, slow relaying of stimuli
- b) Myelinated nerve fibers with saltatory, fast relaying of stimuli

MOTOR (myelinated)

NERVE FIBERS (skeletal muscles)


Up to 10 Hz	Single twitches
10 – 20 Hz	Fibrillations
20 – 80 Hz	Tetanic contractions
above 90 Hz	Muscle waves

AUTONOMIC (unmyelinated)

NERVE FIBERS (smooth musculature, organs)

Below 10 Hz	Stimulation of the vascular muscles
10 – 20 Hz	Stimulation of the intestinal muscles*

* smooth musculature (muscles of the internal organs, deep lying – involuntarily contracting muscles)

INTRODUCTION

Treatment theories

1. Endorphin theory

At very low frequencies (< 10Hz), stimulation current predominantly causes hormonal manipulation of vegetative centers in the central nervous system. Pain-transmitting fibers are stimulated, triggering central pain-inhibiting mechanisms and opioid release in the brain stem or diencephalon, and the body's own endorphins and morphine-like substances are released. These work like painkillers!

2. Gate control theory

Stimulation frequencies > 10Hz have their effect mainly in the "dorsal horns" of the spinal cord. There, pain sensations are prevented from penetrating to consciousness. Frequencies from 85 – 100 Hz to 200 – 250 Hz can trigger pain-inhibiting mechanisms in the free nerve endings and sympathetic nerve fibers.

This is the most widespread theory.

General pain analysis

1. Intense, sharp or piercing pain

Use frequencies < 50 Hz to < 20 Hz



Select as high a pulse width as possible > 150µsec

2. Dull, throbbing pain

Use frequencies > 50 Hz (typically)




Select as wide a pulse width as possible (the pain or throbbing must not be aggravated), that usually means a small pulse width of < 120µS

Avoid stimulating muscles at the same time, otherwise this can cause aching muscles!

How TENS works

Tens electrical pulses override the pulses from the source of pain.

POSSIBLE TENS APPLICATIONS

Scientifically founded, TENS proves effective in treating the following pain conditions:

1. Acute pain conditions:

- Postoperative pain (**P08**)
- Fracture pains (**P09**)
- Childbirth pains (**PC1 or PC2 - 120Hz, $\mu S > 100$**)

2. Chronic pain conditions:

- Pain in the neck, back and small of the back (**P02, P09**)
- Painful stiffness in the shoulders (**P09**)
- Rheumatic complaints (**P08**)
- Ischialgia (**P02, P08**)
- Arthrosis joint pains (**P02, P08, PC1 or PC2**)
- Tendovaginitis (**P09, PC1 or PC2**)
- Pain from Bekhterev's arthritis (**P08, P09**)
- Extremity pains (**P08, PC1 or PC2**)

*** PC1 or PC2:**

120 Hz, select pulse width (μS) according to the part of the body:

- | | |
|-------------------------|-------------------|
| - Face | 50 μS |
| - Shoulder/upper arm | 150 μS |
| - Lower arm | 200 μS |
| - Thorax/thoracic spine | 250 μS |
| - Stomach/lumbar spine | 300 μS |
| - Thighs | 350 μS |
| - Lower leg/foot | 400 μS |

Cervical spine syndrome / myalgia (neck pain)

(P08, P09, PC1 or PC2 - 120Hz, 50 μS)

Can, among other things, be a sign of wear and tear in the cervical spine area, impaired blood circulation or also psychological factors such as stress. Relief can be achieved by relaxing the cervical vertebra zone of the upper trapezoid muscle and the muscle bands of the spine.

Ischialgia

(P08, PC1 or PC2 - 120Hz, 300 μS)

Very good results are achieved when treating back pain and pain in the small of the back. One in two patients has experienced such symptoms at some point. The electrodes can be positioned in the back area or also on the lower extremities.

D

US

F

I

E

Shoulder pain

(P10, PC1 or PC2 – 120Hz, 150µS)

One-sided overexertion or carrying heavy loads, infections and psychological strains.

Hip joint pain

(P08, PC1 or PC2 – 150Hz, 350µS)

Therapy goals are pretty much identical for rheumatoid arthritis, arthrosis, post-traumatic conditions and artificial limbs (only after consulting a doctor): relieving pain and reducing inflammation and muscle tension. Arthrosis, arthritis (knee joint pain). Possible causes for knee joint pain include, for example, metabolic disorders, overweight, mechanical strains and the normal ageing process of the cartilage tissue.

Foot joint pain


(P08, PC1 or PC2 – 120Hz, 400µS)

In general, TENS can be used to treat both chronic and acute pain, e.g. in the case of sprains and dislocations that are usually accompanied by bruising. However, a careful examination is necessary first, to establish whether the joint may be damaged.

Muscle tension

(P01, P03, P09)

caused by, for example, jogging or after long periods of driving or sitting at the computer, can be most effectively treated using burst mode (pulse groups are "fired off").

Thoracic syndrome

(P08, PC1 or PC2 – 120Hz, 250µS)

Thoracic syndrome is not usually a serious condition, but rather a painful muscle cramp, which can be exacerbated by poor posture, activities carried out in a sitting position and depression.

Elbow pain

(P01, P08, PC1 or PC2 – 120Hz, 150µS)

Elbow pain may be caused by one-sided exertion or overexertion of the elbow joint, for example when playing golf or tennis. The problem is an inflammation of the tendons of the extensor muscle. Apart from resting the arm, TENS treatment is an ideal aid in achieving rapid pain relief.

Vein problems / lymph drainage




(P01, P02)

Here, both motor (fast) and autonomic (slow) muscle fibers should be stimulated – using a frequency of approx. 150 Hz reaches both types.

Migraine

(P02, P09)

As the causes of headache vary greatly depending on the individual, various combinations of stimulation points may be tried out. With a little patience, you will soon find the optimum points of stimulation. You should only consider treatment of the head region after consulting your doctor.

POSSIBLE TENS APPLICATIONS

Allergies

(P02, PC1 or PC2 - 120Hz, 150µS)

Position the electrodes on the following master points:

Electrode 1:


At the end of the elbow fold or between thumb and index finger

Electrode 2:

4 patient's finger-breadths below the wrist fold on the outer side of the arm

It is imperative to consult a doctor BEFOREHAND in the following cases:

- Prolonged pain symptoms despite therapy
- When using strong painkillers or local anesthesia
- Infectious diseases
- Circulation problems, thrombosis
- Impaired sensibility (numbness)
- Use with babies and infants
- Pregnancy
- Psychosis
- Bleeding tendency
- Cancer
- Extreme sensitivity to or fear of electricity
- Pacemaker
- Artificial joints or metal plates

TENS THERAPIE RECOMMENDATIONS

Therapy goal

Skin stimulation
Analgesia (pain relief)
Detoning (relaxing, loosening)
Trophic (tissue nutrition)

Indications:

Myalgia (muscle pains)
Loosening up muscles
Muscle pump (set intensity for both channels alternately high and low)
Anal incontinence (sphincter muscle weakness)
Pelvic floor weakness
Sympathicus suppression (suppression of tension in the sympathetic trunk)
Promoting resorption
Joint pains
Mackenzie zones (muscle stimulation points)
Head's zones (sensitive or painful skin zones)

For sensitive persons
For thick muscle groups or athletes

Frequency

approx. 50 Hz
approx. 50–180 Hz
approx. 100–120 Hz
approx. 100–120 Hz




approx. 100–120 Hz
approx. 10 Hz
approx. 50 Hz
approx. 1–2 Hz
approx. 25 Hz
approx. 120 Hz
approx. 100 Hz
approx. 50–100 Hz
approx. 120 Hz
approx. 120 Hz

pulse width < 150 μ s
pulse width > 150 μ s

Alternative

approx. 1–4 Hz
approx. 2–4 Hz
approx. 2–8 Hz

approx. 2–4 Hz

HOW DOES EMS (ELECTRICAL MUSCLE STIMULATION) WORK?

EMS stands for **E**lectronic **M**uscle **S**timulation.

Muscle contractions are triggered by electrodes positioned on the body. The successive contracting and relaxing of the muscles strengthens them and at the same time helps integrate the muscles more quickly into general movement patterns, promoting and stimulating both muscle activity and arterial circulation. The **Promed EMT-4** is ideally suited to preventing and delaying disuse atrophy in injured parts of the body, e.g. when a part of the body is immobilized for a long period of time due to a plaster cast, splint or extension bandage.

A further area of application for the **Promed EMT-4** is the treatment of various mobility disorders, for example involving the hip, ankle and elbow joints. The stimulator can be used on the corresponding parts of the body to improve coordination in the case of mobility dis-

orders. This acts on both the **nerve** and **muscle tissue**. By loosening the muscles, pain can be relieved or even eliminated. **The stimulation frequencies are lower than with the TENS method.**

During the entire cycle time, the muscles being treated can be

- **tensed**
- **relaxed**
- **moved**
- **kneaded**
- **massaged**

thereby achieving a clear reduction in local pain.

Both the TENS and EMS methods can activate the muscles and thus have a muscle-building effect.

Low frequencies (1-10Hz)

with long pulse times and single muscle contractions provide:

- detoxification
- relaxation
- improved circulation in the muscles
- carrying away metabolic end-products (lymph drainage)

Medium frequencies (20-50Hz)

The quick consecutive contractions (fibrillation) provide:

- increased strain on the muscles
- increased muscle building
- increased muscle fiber size
- improved oxygen uptake
- increased energy consumption for the user

D

US

F


I

E

POSSIBLE EMS APPLICATIONS:

Applications

- Regenerating parts of the body that have been at rest for long periods (plaster cast)
- Treating other mobility disorders, e.g. with the hip, ankle and elbow joints
- Supporting training performance and muscle building
- Warming up the muscles
- Muscle stamina training
- Muscle relaxation
- Muscle detox
- Improving circulation and metabolism
- Rehabilitation after injuries
- Pain therapy
- Stimulating arterial circulation
- Promoting muscle activity, especially in older people
- Zur Koordinationsverbesserung gestörter Bewegungsabläufe
- Improving coordination in the case of mobility disorders
- Preventing thrombosis
- Treating scoliosis
- Cellulitis/toning: this shows how versatile **Promed** electrotherapy can be. With regular use, combined with an exercise and nutritional program, you can achieve visible success. All TENS/EMS pulse current types are suitable here. Several treatments a day are recommended.

MASTER POINTS AS A SOURCE FOR IMPROVING PERFORMANCE




These points, stemming from acupuncture, can be used at any time as support or can even be chosen as main points. The master points are to be found on the lower arm, chest and leg.

Lower arm master points (Pe 6, between the tendons of the radial flexor muscles of the wrist on the inside of the lower arm, 2 Cun proximal to the wrist and Lu 9, on the radial side of the fold between hand and wrist, lateral to the Arteria Radialis) have a relaxing effect on the cardiovascular system, and improve circulation and blood flow properties.

Chest master points (Lu 1, one Cun below the collar bone, six Cun lateral to the center line between the first and second ribs) increase breathing capacity.

Leg master points (Ma 36, one finger's breadth lateral to the lower edge of the tuberositas tibia/tibial plateau, three Cun below the knee fold) are said to be the best general tonifying points.

The master points can be used individually or in combination. They are present on both the left and right sides of the body.


CONTRAINDICATIONS

Before using the unit, always be sure to consult a doctor, especially in the following cases!

- Patients with implanted electronic devices (e.g. pacemaker)
- Patients with metal implants
- During pregnancy
- With children
- With heart problems, especially arrhythmia
- Before any positioning of the electrodes that transmits electricity to the carotid sinus nerve (neck)
- Before any positioning of the electrodes that transmits electricity transcerebrally (through the head)
- With non-diagnosed pain symptoms
- Treatment on the eyelids




WARNINGS

- **Promed EMT-4** treatment must not be carried out at the same time as another potentially dangerous activity, e.g. driving an automobile, operating a machine etc.
- It must not be used during bathing or showering.
- The **Promed EMT-4** must be kept safely out of the reach of children.
- Do not misuse the unit or use it for purposes other than the designated applications.
- Protect the **Promed EMT-4** from moisture.
- The **Promed EMT-4** and its accessories must only be used as detailed in the operating instructions.
- You must follow the instructions supplied with the electrodes. We advise against the use of electrode sizes other than those recommended by the manufacturer.




PRECAUTIONARY MEASURES

- **Promed EMT-4** must not be used in the vicinity of microwave devices, short wave devices, high-frequency surgical equipment or in airplanes.
- We also warn specifically of the following dangers:
 - Simultaneously connecting the patient to a high-frequency surgical device can result in burns under the stimulation current electrodes.
 - Using the unit near (e.g. 1 m) a short wave or microwave therapy device can cause fluctuations in the output of the stimulation current.
- Should treatment with the **Promed EMT-4** be ineffective or cause discomfort, discontinue use and consult a doctor.
- Make sure you do not use the unit if it or the cables are damaged or faulty. If a fault occurs, please send in the unit or cables via your specialist dealer to:
Promed GmbH
Service Department, Lindenweg 11, D-82490 Farchant.
- Long-term use can in isolated cases lead to skin irritation in the area where the electrodes were placed.
- The effectiveness of the treatment depends greatly upon the patient being treated by someone who is qualified in dealing with patients suffering from pain.


POSITIONING ELECTRODES

DISPLAY GUIDE


When the unit is switched on, all readouts are active for a short time before the unit switches to standby (wait) mode.

PROGRAM SELECTION...

Preset programs P-01 to P-11 and memory function

1. Press the "ON/OFF" button to switch on the unit.
2. Keep pressing the "Program" button until one of the preset programs (P-01 to P-11) is shown on the display.
3. Press the "Set" button. The time will flash on the display, " ' " for minutes and " ` " for seconds. Use the "+" or "-" buttons to enter the desired treatment duration (20 minutes to 2 hours) in steps of ten minutes.
4. Press the "Select" button 2 times. The intensity display "CH1" will start to flash. Use the "+" or "-" buttons to enter the desired value (from 1 to 50). Press the "Select" button once more. The intensity display "CH2" will start to flash. Once again, use the "+" or "-" buttons to set the desired value.
5. To begin treatment, press the "Select" button once more. The clock symbol at the top right of the display will start to flash.

AND STARTING THE TREATMENT

MEMORY FUNCTION

6. If you want to change the intensity during treatment, press the "Select" button again and adjust the intensity as required using the "+" or "-" buttons.




Calling up and resetting the memory

1. Press the "ON/OFF" button to switch on the unit.
2. Press and hold down the "Set" button for more than 3 secs. to switch to the memory function. The total treatment time will now be displayed for 30 secs.

After this, the unit will switch back to standby mode, or you can press any button to return to setting mode.

3. When in memory mode, hold down the "Select" button for 3 secs. to reset the treatment time to zero.

If no buttons are pressed for at least 2 minutes before confirming a function, the unit will automatically return to time display mode. To switch off the unit, press the "ON/OFF" button once more.

PROGRAM SELECTION AND STARTING THE TREATMENT

Programmable programs P-C1 and P-C2


1. Press the "ON/OFF" button to switch on the unit.
2. Keep pressing the "Program" button until the programmable programs (P-C1 or P-C2) are shown on the display.
3. Press the "Set" button once. In the display, the frequency (Hz) and "SET" will start to flash. Use the "+" or "-" buttons to enter the desired value (1 Hz to 200 Hz).
4. Press the "Set" button. In the display, the pulse width (μ S) and "SET" will start to flash. Use the "+" or "-" buttons to enter the desired value (from 30 μ S to 400 μ S) in steps of ten.
5. Press the "Set" button once more. The time will flash on the display, " ' " for minutes and " ` " for seconds. Use the "+" or "-" buttons to enter the desired treatment duration (1 minute to 24 hours).
6. Press the "Set" button to confirm and then press the "Select" button 2 times. The intensity display "CH1" will start to flash. Use the "+" or "-" buttons to enter the desired value (from 1 to 50). Press the "Select" button once more. The intensity display "CH2" will start to flash. Once again, use the "+" or "-" buttons to set the desired value. Fast adjustment protection will be activated and stops at intensity level 15 when the intensity increases rapidly.
7. To begin treatment, press the "Select" button once more. The clock symbol in the top right of the display will start to flash.
8. If you want to change the intensity during treatment, press the "Select" button again and adjust the intensity as required using the "+" or "-" buttons.

Pause mode

1. To temporarily pause treatment, press the "Program" button. In the display, the pause symbol will start to flash. To continue, press the "Program" button once more.
2. To stop treatment ahead of time, press the "ON/OFF" button 2 times. The unit will revert to standby mode and display the time.
3. Once the treatment period is completed (countdown), the unit will automatically switch to standby mode.


Setting the time

1. In standby mode, the unit display the time. If it is showing the treatment settings, press the "ON/OFF" button.
2. To set the clock, now press the "Set" button. The hours will start to flash. Use the "+" or "-" buttons to set the desired hour. Press the "Set" button once more. The minutes will start to flash. Use the "+" or "-" buttons to set the minutes. Press the "Program" button to confirm.




PROMED EMT-4 PROGRAM TABLE

Programm	Pulse 1		Pulse 2		TENS/EMS	Quantity per cycle	Which program for which purpose?
	Frequency (Hz)	Pulse width (µS)	Frequency (Hz)	Pulse width (µS)			
P01	15	200	Burst frequency = 2 Hz		TENS	Burst, 54 times, 56 secs. break 5 secs. per cycle	Extremities (legs, arms) and back, warming up, loosening up, relaxing
P02	150	30	Burst frequency = 2 Hz		TENS	Burst, 54 times, 29 secs. break 5 secs. per cycle	Face, headaches, promotes concentration
P03	5	100	Burst frequency = 2 Hz		EMS	Burst, 3 times, 5,8 secs. break 11,6 secs. per cycle	Loosening up, stamina, aerobic metabolism
P04	75	100	60	300	EMS	Modulation, 5 secs. break 10 secs. per cycle	Muscle building, hypertrophy, anaerobic metabolism
P05	50	175	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Strength endurance, anaerobic metabolism, rehabilitation
P06	110	175	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Arms, maximum strength, intramuscular coordination
P07	100	125	100	225	EMS	Modulation, 5 secs. break 10 secs. per cycle	Arms, legs, maximum strength, intramuscular coordination
P08	130	100	-	-	TENS	Conventional, 31 times break 5 secs. per cycle	Extremities, back, shoulders, pain relief
P09	35	30	Burst frequency = 2 Hz		TENS	Burst, 40 times, 29 secs. break 5 secs. per cycle	Neck, face, shoulders, pain relief, loosening up
P10	40	100	-	-	TENS	Conventional, 62 secs. break 5 secs. per cycle	Shoulder pain, knee complaints
P11	50	250	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Breast toning, pectoralis training, rehabilitation
PC1	1-200	30-400	-	-	TENS	-	Individual treatment, free to be set and stored by patient or doctor
PC2	1-200	30-400	-	-	TENS	-	Individual treatment, free to be set and stored by patient or doctor

TECHNICAL SPECIFICATIONS


- Power supply: 4 x batteries 1,5 V, AAA type
(total 6,0 V) or 4 x rechargeable batteries 1,2 V AAA type
- 2 channels: CH 1 and CH 2
- Current type: biphasic square pulse
- Intensity: 0 mA – 50 mA, at 500 mW
(each channel is separately adjustable, 50 stages)
- Frequency range: 1 Hz to 200 Hz
- Pulse width range: 30µsec to 400 µsec
- treatment time:
 - a) Preset programs: (P-01 to P-11):
20 mins to 2 hours
 - b) Programmable programs: (P- C 1 and P- C 2):
1 min to 24 hours
- Environmental conditions:
 - Temperature: 5° to 40° C
 - Relative humidity: 30 to 80 %
- Storage conditions:
 - Temperature: - 10° to + 55° C
 - Relative humidity: 10 to 90 %
- Dimensions: 122 x 58 x 18 mm (L x W x H)

FURTHER INFORMATION

Certification No. CE 0123

Classification label

V
mA
mW
Hz
 μ s
k Ω m
mm
g
BF
P


Volts
Milliamps
Milliwatts
Hertz
Micro second
Kilo-Ohms
Millimetres
Grammes
Body floating
Output rating

General warning.
According to DIN VDE 750/219, the unit must be marked with this symbol. It means that the user must follow these operating instructions and thus use the unit for the purpose intended.

Direct Current

TROUBLESHOOTING

SERVICE AND MAINTENANCE

If the Promed EMT-4 is not functioning correctly

1. Make sure that the batteries are correctly inserted, or replace the batteries. When replacing the batteries, make sure that the polarity is correct.
2. If the display is showing all functions but you can not feel any stimulation, check that the electrode cables are properly connected and the electrodes are correctly positioned. If the unit seems to be working, but you feel no stimulation, it is possible that the electrode cables or the electrodes have to be replaced.
3. If the unit seems to be working, but you still feel no stimulation, increase the current intensity as described above.
4. If you still can't feel any stimulation, the electrodes or electrode cables may have to be replaced.

Should you encounter further problems, please contact your dealer or return the unit to your dealer. Do not attempt to repair a faulty unit yourself.




For repair or maintenance send the **Promed EMT-4** to the following address via your specialist dealer:


Promed GmbH

Service Department

Lindenweg 11, D-82490 Farchant

Tel. +49 (0) 8821 96 21-19, Fax +49 (0) 8821 96 21 21

CLINICAL RECOMMENDATIONS

Clinical recommendations for Promed EMT-4 settings

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









US | Guarantee conditions

The manufacturer provides a 24 month guarantee for the **Promed EMT-4** stimulator, valid as of the date of purchase shown on the receipt or invoice. This guarantee does not apply to the accessories (electrodes, leads etc.) The guarantee does not cover: damage caused by force, external factors or influences, damage caused by incorrect use, any repairs, maintenance or modifications carried out by any unauthorized third party and damage resulting from not reading the instructions.

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